


H HOTPLATE/STIRRERS

Multi Position

Multi-Hotplate Stirrers, 4 Or 6 Places

Aluminium Top Plate • Independent Heating & Stirring Control • up to 300°C/1600rpm

Excellent for Multi-use & a Wide Various Sample Testing • Systemized 4 & 6 Hotplate-Stirrer • Independent Heating & Stirring Control • Maximum Temperature of 300°C • Aluminium tops heat up quickly & uniformly. Great for low temperature applications • Units are built from durable cast aluminium cases & painted with chemically resistant paint.



MUHS-6

Model: MUHS-4, 4 Positions

Multi-position hotplate stirrer with four positions (in one row), Same as model MUHS-6, except dimensions which

are: 60.5x14.5x22cm, weight 8 kg, power 2000 watt-4x500 watt.



MUHS-2x2



MUHS-3x2

Model: MUHS-6, 6 Positions

- Multi position hotplate stirrer with six positions.
- Six individually controlled stirring and heating positions.
- Each position can be used to stir only, heat only or heat and stir at the same time.
- Ideal for a narrow bench-top or for applications needing support rack system.
- The top plates are chemical resistant and heat quickly.
- Stirrers are equipped with strong magnets & high torque motors.

Temperature range: Up to 300°C

Speed control: Up to 1600RPM

Dimensions: Top plate 11.5x11.5cm,
Overall W81xH14.5xD22cm
(D includes support rod clamp)

Power: 230V (or 115V), 50/60Hz, 3000W-6x500W

Net weight: 12Kg

300°C

Options: Model: MUHP-4

Similar to the above described model: MUHS-4, but without stirrers.

Model: MUHP-6

Multi position hotplate with six positions.

Model: MUHS-2x2

- Multi position hotplate stirrer with four positions (in 2 rows 2+2).
- 4 individually stirring and heating positions.
- 4 chemical resistant top-plates, 11.5x11.5cm.
- The top plates are chemical resistant and heat quickly.
- Stirrers are equipped with strong magnets and high torque motors.

Temperature range: Up to 300°C

Speed control: Up to 1600RPM

Dimensions: Top plate 11.5x11.5cm
Overall W30xH14.5xD31cm
(D includes support rod clamp)
Placed in two rows 2+2

Power: 230V (or 115V), 50/60Hz, 2000W-4x500W

Net weight: 8Kg

Options: Model MUHP-2x2, Similar to the above described model: MUHS-2x2, but without stirrers

Model: MUHS-3x2

Multi-position hotplate stirrer with six positions (in 2 rows 3+3)

Same as model MUHS-2x2, except dimensions which are:

45x14.5x31cm, weight 8 kg, power 2000 watt-4x500 watt.

- 6 individually stirring and heating positions.
- 6 chemical resistant top-plates, 11.5x11.5cm.
- The top plates are chemical resistant and heat quickly.
- Stirrers are equipped with strong magnets and high torque motors.

Options: Model MUHP-3x2, Similar to the above described model: MUHS-3x2, but without stirrers