

CENTRIFUGES 300ml Capacity



SCEN-207, Beginner's Small Universal Centrifuge

Features:

- Analog speed tachometer.
- Optional rotor selections for various tube sizes.
- Regressive timer: 30min, or continuous operation.
- Angle rotor, stepless speed setting.
- Auto accelerating and decelerating functions.
- Auto-current restraining function.
- Power pilot light and fuse safety device.
- Manual cover lock.
- Speed: 5,500rpm; Timer: 0~30 or "ON".
- Max. capacity: 300 ml.
- Weight: 9kg.
- Size: W29xD32xH25cm.
- Auto power off when open and auto-balance device.
- ABS hard shell and see-through window.
- Auto-cooling air circulation.



Model	SCEN-207
Max. Speed (rpm)	5,500
Max. RCF (g)	3,280
Max. Capacity (ml)	300
Timer (min)	0~30 or "ON"
Weight (kg)	9
W x D x H Size (cm)	29 x 32 x 25

Rotors for SCEN-207



RA-1508
15ml x 8pcs
Max. 5,500rpm



RA-1512S
15ml x 12pcs
Max. 5,500rpm



RA-5004
50ml x 4pcs
Max. 5,500rpm



REA-1520
1.5ml x 20pcs
Max. 4,000/5,500rpm



RA-5006
50ml x 6pcs
Max. 5,500rpm

	Max. Speed (rpm)	Max. RCF (g)	Capacity (ml x pcs)	Tube size (Ø x L) mm
RA-1508	0~5,500	3,140	15 x 8	15~16.7x90~122
RA-1512S		3,280	15 x 12	15~16.7x90~122
RA-5004		3,240	50 x 4	27.5~29.5x90~117
RA-5006		3,240	50 x 6	27.5~29.5x90~117
REA-1520		2,600	1.5 x 20	10~11x36~56



SCEN-206/HSCEN-204

Model	SCEN-206	HSCEN-204
Max. Speed (rpm)	500~6,000	1,000~10,000
Max. RCF (g)	3,860	10,390
Max. Capacity (ml)	300	
Timer (min:sec)	99:59	
Weight (kg)	10	
W x D x H Size (cm)	29 x 32 x 25	

SCEN-206/HSCEN-204, Digital Universal Centrifuges

Features:

- LED speed and time display.
- Microprocessor control.
- 10 memory programs for preset parameters.
- Auto accelerating and decelerating functions.
- Auto cover lock with motor running.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- Speed: 500~6000rpm ;
Timer: 99min 59 sec.
- Max. Capacity: 300ml.
- Weight: 10kg.
- Size: W29xD32xH25cm.
- ABS hard shell and see-through window.
- Auto-cooling air circulation.
- Optional rotor selections for various tube size.



Rotors for SCEN-206



RA-1508
15ml x8pcs
Max. 6,000rpm



RA-1512S
15ml x12pcs
Max. 6,000rpm



RA-5004
50ml x4pcs
Max. 6,000rpm



RA-5006
50ml x6pcs
Max. 6,000rpm



REA-1520
1.5ml x20pcs
Max. 6,000rpm

Rotors for HSCEN-204 (High-Speed)



RA-1508
15ml x8pcs
Max. 10,000rpm



RA-1512S
15ml x12pcs
Max. 9,000rpm



RA-5004
50ml x 4pcs
Max. 8,000rpm



RA-5006
50ml x 6pcs
Max. 7,000rpm

Rotor	HSCEN-204		SCEN-206		Capacity (ml x pcs)	Tube size (Ø x L) mm
	Max. Speed (rpm)	Max. RCF (g)	Max. Speed (rpm)	Max. RCF (g)		
RA-1508	10,000	10,390	6,000	3,740	15 x 8	15~16.7x90~122
RA-1512S	9,000	8,780	6,000	3,900	15 x 12	15~16.7x90~122
RA-5004	8,000	6,860	6,000	3,860	50 x 4	27.5~29.5x90~117
RA-5006	7,000	5,250	6,000	3,860	50 x 6	27.5~29.5x90~117
REA-1520	–	–	6,000	3,090	1.5 x 20	10~11x36~56