



BSCEN-600

BSCEN-600/601/602, Centrifugal Separator (Tripedal Type)

Features:

- For solid and liquid separation.
- The dipping part is in stainless steel SUS304.
- For materials separation in large volume.
- Equipped with 3-point shake-proof device.

Teflon coating for getting acid resistance can be requested to do with an extra charge.


BSCEN-600




**BSCEN-600T
Teflon Coated**



Filter bag



**Metal Basket
(Small)
Ø15cm x H8cm
(Capacity 0.8L)**



**Metal Basket
(Medium)
Ø24cm x H14cm
(Capacity 3L)**



**Metal Basket
(Large)
Ø30cm x H16cm
(Capacity 4.5L)**

Model	Max. Speed (rpm)	Max. RCF (g)	Basket Size (cm)	Capacity (liter)	Weight (kg)	W x D x H Size (cm)
BSCEN-600	0~3,000	1,509	Ø15 x H8	0.8	74	56 x 56 x 60
BSCEN-601			Ø24 x H14	3.0		
BSCEN-602			Ø30 x H16	4.5		