



OLCEN-501/502

OLCEN-501, Oil Test Typed Centrifuges



Model	OLCEN-501	OLCEN-502
Max. Speed (rpm)	500~2,000	
Max. RCF (g)	1,070	1,140
Max. Capacity (ml)	4	8
Timer (min:sec)	99:59	
Weight (kg)	160	
W x D x H Size (cm)	69 x 77 x 104	

Tubes not included


AC110/220V, 50/60Hz

Features:

- LED speed, centrifugal force and time display
- 10 memory programs for presetting parameters.
- Auto rotor, stepless speed setting
- See-through window.


This machine can be applied to test the water content and sediment in the oil or the insolubility insolubility ingredient in the lubricant.

Rotor Included



Oil Test Bottles

OLCEN-501



RS-100

Tube x 4
GT-1/GT-2
GT-4/GT-5

OLCEN-501/502



RS-104

GT-3
(Sole GT-3 for RS-104)

OLCEN-502



RS-108

Tube x 8
GT-1/GT-2
GT-4/GT-5

Rotor for OLCEN-501/502

Rotor	Max. Speed (rpm)	Max. RCF (g)	Capacity (ml x pcs)	Tube size (Ø x L) mm
RS-100	2,000	1,070	Refer to the above diagramatic explanation	
RS-104		850		
RS-108		1,140		