



HSCEN-1001

Features:

- LCD, speed and RCF setting also display.
- 10 memory programs for presetting parameters.
- Auto accelerating and decelerating functions.
- Auto power off for rotor imbalance, overload, breakdown indication and safety protection device.

HSCEN-1001, High Speed Centrifuge



Model	HSCEN-1001
Max. Speed (rpm)	1,000~10,000
Max. RCF (g)	10,950
Capacity (ml x pcs)	400
Timer (min:sec)	99:59
Weight (kg)	35
W x D x H size (cm)	36 x 48 x 34

Tubes not included

AC110/220V, 50/60Hz

Rotors for HSCEN-1001



RA-1508
15ml x 8pcs
Max. 6,000rpm



RA-1512H
15ml x 12pcs
Max. 6,000rpm



RA-5004
50ml x 4pcs
Max. 6,000rpm



RA-5006
50ml x 6pcs
Max. 6,000rpm



RA-5008
50ml x 8pcs
Max. 6,000rpm

Rotor	Capacity (ml x pcs)	Tube size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (g)
RA-1508	15x8	15~16.7x90~122	10,000	10,950
RA-1512H	15x12	15~16.7x90~122	9,000	10,130
RA-5004	50x4	27.5~29.5x94~130	8,000	7,080
RA-5006	50x6	27.5~29.5x90~130	7,000	5,400
RA-5008	50x8	27.5~29.5x90~117	6,500	5,280